



AirFit™ F30i

Fitting guide for AirFit F30i



1 With one hand holding the bottom of the cushion and the other holding the conduit frame, position the mask in front of the face.



2 Place the cushion under the nose and against the face. With the ResMed logo and headgear facing up, pull the headgear over the head.



3 Bring the lower headgear straps under the ears and attach the magnetic clips to the frame connectors.



4 Undo the fastening tabs on the upper headgear straps and pull evenly.



5 Repeat with the lower headgear straps.



6 Squeeze the side buttons on the elbow and detach from the frame. Connect the air tubing from the device to the elbow, then reattach it to the frame.



For people with long hair, we recommend feeding the hair through the back of the headgear, as pictured. If tied in a ponytail, it should sit high enough to fit through the back of the headgear.

Refer to the user guide for further instructions and complete labelling information including indications, contraindications, warnings and precautions.

¹ ResMed onsite fitting study of 73 current full face CPAP mask users and 14 sleep professionals from Nov. 5-14, 2018 in multiple US locations who trialled AirFit F30i and a competitor full face mask.

² Use of masks with magnetic components is contraindicated in patients with the following pre-existing conditions: a metallic hemostatic clip implanted in the head to repair an aneurysm and/or metallic splinters in one or both eyes.

³ ResMed guided external clinical study of 48 evaluable current full face CPAP mask users from 24/10/2018 to 20/11/2018 in Sydney, AU and San Diego, CA, who trialled AirFit F30i and competitor full face mask at home. Data on file; ID A4484728.

© 2020 ResMed Ltd. 10111833 2020-02

ResMed.com/AirFitF30i



AirFit™ F30i



The new face of bedtime freedom

Give your customers a mask that lets them do more in bed when awake or asleep. ResMed AirFit™ F30i is a tube-up UltraCompact full face mask designed with a reliable seal, under-the-nose cushion and the freedom to sleep in many different positions. In a ResMed study, users rated AirFit F30i 4/5 stars for overall mask comfort and seal.¹



Awaken your best.


ResMed.com/AirFitF30i

AirFit F30i highlights

Quick-release elbow is designed for easy connection and disconnection of the mask from tubing, while the 360° elbow rotation allows for freedom of movement.

Flexible conduit frame allows therapy air to flow through the elbow and into the frame and cushion.

UltraCompact full face cushion rests under the nose to deliver all-round stability¹, visual freedom, and minimal facial contact.




Top-of-head tube design gives users the freedom to sleep the way they want.

Magnetic clips make it quick and easy to connect the headgear to the frame.²

Optional soft wraps can help add softness and comfort at the cheeks.

Made for a natural sleep experience

ResMed's freedom-focused AirFit F30i allows customers who toss and turn at night to move comfortably in their sleep.³ Compared to a traditional full face mask, the tube-up design and 360° rotation of AirFit F30i offer users a more natural sleep experience. If your customer is an active sleeper, give them a full face mask that lets them sleep in many different positions.



Product codes

AirFit F30i: S/SML	ANZ 63336	AirFit F30i: M/STD	63338
AirFit F30i: S/STD	63337	AirFit F30i: W/STD	63339
		AirFit F30iL M/LGE	63340

Awaken your best.

The UltraCompact, under-the-nose cushion of AirFit F30i full face mask means there's no discomfort or red marks on the bridge of the user's nose.



AirFit™ F30i
Full face mask