



Live life in high Res

There's a full face
mask for every
sleeper



Drift to sleep wearing our most comfortable gear

Let's find a mask that fits your needs

At Resmed, we believe a mask should seal securely and comfortably—no matter how you sleep or breathe. That's why each of our masks is designed with different features to meet a variety of needs. After all, every person and every face is unique. With the number of options we have available, we're confident you'll discover one that feels tailor-made.



What is a full face mask?

A full face mask typically has a cushion that covers both your mouth and nose to create a seal over both airways.

If you sleep and breathe with your mouth open or experience nasal congestion, you'll be recommended a full face mask. They are also more suitable for therapy at higher pressures as they are more tolerable than other mask types.



We're confident you'll love our full face masks — and if you're not satisfied, we've got you covered with a 30 day money-back guarantee or a 14 day cushion replacement guarantee. Just check the icons to see which masks qualify.*

*Terms and conditions apply. See in store for details.

Rest assured, we'll help you find the right mask

You're sorted with our mask categories

Our full face mask range is continually growing and evolving to get better. With our mask categories, finding the right mask has never been easier. These categories make it easy for you to understand the key benefits of each mask, so you can take one home with confidence. Each category offers something a little different, so have a look to see which one suits you best.

Minimalist masks

Small. Light. Simple.

Freedom masks

Sleep and move comfortably.

Versatile fit masks


Designs that fit nearly every face.^{1,2,3,4}

Speak to your Resmed Authorised Dealer to find the right mask that will help you adjust to therapy, and receive optimum results night after night.



STANDARD
EASGEAR

ResMed



Minimalist

AirFit F40 is Resmed's smallest UltraCompact full face mask with a flexible design that allows you to sleep in your favourite positions.

AirFit F40 is available in SW/STD, M/SML, M/STD and L/STD sizes.

- 1 UltraCompact AdaptiSeal™ cushion made of a soft, total silicone rests securely under the nose and adapts to your face.
- 2 Flexible curved mask frame sits below the eyes, designed for comfort and to reduce facial marks.
- 3 QuietAir™ diffused venting gently disperses exhaled air for a more peaceful bedroom environment.
- 4 Quick-release elbow with short tube makes it easy to connect and disconnect to therapy while also reducing tube drag.
- 5 Magnetic clips easily snap together for fast fitting.¹



References

¹The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labelling information, including magnet contraindications and warnings.

AirFit™ F30i

30 DAY MONEY BACK GUARANTEE



Did you know?

This mask is also available without magnetic clips.

Freedom

AirFit F30i is an UltraCompact full face mask that helps keep air tubing up, out of sight and out of the way. It's designed to let you move at night and sleep in many positions while achieving a reliable seal.¹ It's also ideal if you're looking for visual freedom in bed.

AirFit F30i is available in S/SML, S/STD, M/STD, W/STD and M/LGE. Non-magnetic option available.

- 1 Flexible conduit frame allows therapy air to flow through the elbow and into the frame and cushion.
- 2 UltraCompact cushion sits under the nose, eliminating red marks and discomfort on your nasal bridge.
- 3 Magnetic clips guide the headgear to the frame quickly and easily.^{2,3}



References

¹Resmed onsite fitting study of 73 current full face CPAP mask users and 14 sleep professionals from Nov. 5-14, 2018 in multiple US locations who trialed AirFit F30i and a competitor full face mask.

²The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labelling information, including magnet contraindications and warnings.

³Non-magnetic option available.







Freedom

The AirTouch™ F30i Comfort full face mask features a soft, flexible frame wrapped in breathable, moisture-wicking fabric. It combines innovative AdaptiSeal™ technology within a ComfiSoft™ fabric cushion, surrounding you with comfort at every touchpoint. Designed to elevate your sleep experience, it's made to feel as soft and familiar as slipping into your favourite pyjamas.

AirTouch F30i Comfort is available in SW/SML, SW/STD, M/SML, M/STD, M/LGE, and L/STD sizes.

- 1 UltraCompact ComfiSoft™ fabric cushion with AdaptiSeal™ technology that gently contours to your face
- 2 Flexible fabric-wrapped frame, designed to be breathable and moisture-wicking¹ to keep your skin feeling cool and fresh
- 3 QuietAir™ diffused venting gently disperses exhaled air to help make sleep feel more peaceful for you and your bed partner
- 4 Magnetic Clips guide the headgear to the frame quickly and easily²
- 5 Top-of-head tube connection and quick-release elbow — with no tubing in front and a 360-degree rotation, you can sleep in any position you want



References

¹Resmed internal testing for the moisture-wicking properties of fabric, ID: A5498268


²The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labelling information, including magnet contraindications and warnings.

AirFit™ F20



**MONEY BACK
GUARANTEE**

Did you know?
This mask is also available
without magnetic clips.



Versatile fit

Our AirFit F20 full face masks combine the tried and tested format many users know and trust with our best technologies. With a fit range of 96%^{1-3*}, these masks comfortably fit almost all facial shapes and sizes, so you can rely on these masks for performance.

AirFit F20 is available in Small, Medium and Large sizes. AirFit F20 for Her is available in Small and Medium sizes. Compatible with AirMini. Non-magnetic option available.

- 1 InfinitySeal™ silicone cushion is interchangeable with AirTouch™ F20 UltraSoft™ memory foam cushion.
- 2 QuietAir™ diffuser vent minimises noise and gently disperses exhaled air.
- 3 Magnetic clips guide the headgear to the frame quickly and easily.^{4,5}
- 4 Quick-release elbow makes it simple to disconnect from the mask tubing without removing the mask.



*References on backpage.

AirTouch™ F20



Versatile fit

AirTouch F20 features our UltraSoft memory foam**, which is designed to provide a light and breathable seal. Because the foam is softer than silicone, it may help reduce the risk of red marks and skin irritation. AirTouch F20 also shares the same frame and headgear as AirFit F20, as well as the QuietAir diffuser vent, magnetic clips⁴ and quick-release elbow.

AirTouch F20 is available in Small, Medium and Large sizes.



**14
DAY**

**Cushion
replacement
guarantee**

If you are not completely satisfied with the UltraSoft memory foam cushion, you may exchange the cushions that are a part of the AirTouch F20 Starter Kit for one AirFit F20 InfinitySeal™ silicone cushion within 14 days of purchase.[^]

**Recommended to be replaced monthly. ^T&C's apply. See in store for details.

Become a VIP today...

Join sleepvantage and be treated to a range of exclusive member benefits.*



Members-only warranty**

Additional sleepvantage membership warranty on Resmed products.



Feel supported

Access to our friendly product support line.



Stay informed

Access to a range of educational articles and videos.



Get special offers

Be part of exciting promotions and special offers.

For more information, contact sleepvantage today



sleepvantage.com.au



1300 305 705 (AU) or 0800 737 633 (NZ)



info@sleepvantage.com.au

*sleepvantage membership is available to residents of Australia and New Zealand.

**Terms and conditions apply. See website for more information.

AirFit F20 & AirTouch F20 References

¹ Resmed AirFit F20 internal Australian fitting study of 27 existing CPAP patients, conducted between March–April 2016.

² Resmed AirFit F20 internal USA fitting study of 34 existing CPAP patients, conducted April 2016.

³ Resmed AirFit F20 internal EU and APAC fitting study of 90 existing CPAP patients, conducted June 2016.

⁴ The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labelling information, including magnet contraindications and warnings.

⁵ Non-magnetic option available

ALWAYS READ THE LABEL AND FOLLOW THE INSTRUCTIONS FOR USE.

You should speak to your doctor about your symptoms and whether a CPAP device is suitable for you.



Resmed Pty Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

SCAN
HERE TO
SIGN UP



AirFit, AirMini, AirTouch, ComfiSoft, InfinitySeal, QuietAir, sleepvantage and AdaptiSeal are trademarks and/or registered trademarks of the Resmed family of companies. For patent information, see Resmed.com/ip. Specifications may change without notice.

© 2025 Resmed Pty Ltd. 10111834r5 2025–10