



**ResMed**

# Awaken your best.

Nasal and nasal pillows  
range of masks





## Rest assured, it's easy to get the right mask

Our nasal and nasal pillows mask range will continue to grow and get better, and with our new mask categories, finding the right mask has never been easier. These categories make it easy for you to understand the key benefits of each mask, so you can take one home with confidence. Each category offers something a little different. Have a look to see which one suits you best.

### Minimalist category

Small. Light. Simple.

### Freedom category

Sleep and move comfortably.

### Ultra Soft category

Patented memory foam offers unique comfort.

### Universal Fit category

Classic designs that fit nearly every face.

Speak to your ResMed Authorised Dealer to find the right mask that will help you adjust to therapy, and receive optimum results night after night.

## Drift to sleep wearing our most comfortable gear

Let's find a mask that fits your needs

At ResMed, we believe a mask should seal securely and comfortably – no matter how you sleep or breathe. That's why each of our masks is designed with different features to meet a variety of needs. After all, every person and every face is unique. With the number of options we have available, we're confident you'll discover one that feels tailor-made.

Enjoy a relaxing night's sleep wearing the latest from our nasal and nasal pillows mask range.



NASAL MASK



NASAL PILLOWS MASK

### What is a nasal mask?


Most nasal masks consist of a cushion that covers the nose between the nasal bridge and upper lip. Some nasal masks have cushions that cradle under the nose. Both types are great options if you:

- Breathe through your nose
- Prefer less coverage on your face

### What is a nasal pillows mask?

A nasal pillows mask is a lightweight and compact mask with tubular silicone inserts that fit gently inside your nostrils to deliver therapy air. Nasal pillows masks are a great option for:

- Reading or watching TV in bed
- Reducing facial markings



30  
DAY

MONEY BACK  
GUARANTEE

We're so sure that you will love our AirFit™ nasal and nasal pillows masks that we are offering a 30 day money back guarantee if you are not completely satisfied.\*

\*Terms and conditions apply. See in store for details.

# AirFit™ N30

Nasal cradle mask



## Minimalist

ResMed's lightest mask ever\*, AirFit N30, is designed to make staying on therapy easy. Small and easy to use, it features a curved nasal cradle cushion and adjustable QuickFit™ headgear.

AirFit N30 includes small, medium and small-wide cushions.


- 1 **Adjustable QuickFit headgear** is designed to make fitting easy and adjustment easy.
- 2 **Nasal cradle cushion** curves under the nose for a comfortable and secure fit.
- 3 **QuietAir™ vent technology** diffuses exhaled air flow gently and quietly.



\*As of October 2019. The total weight of AirFit N30 is 45g, inclusive of the headgear but non-inclusive of the packaging.

# AirFit™ P10

Nasal pillows mask




## Minimalist

AirFit P10's sleek styling provides a real sense of facial freedom, and with just three parts, it's designed to be simple to use and maintain. AirFit P10 is our best-selling nasal pillows mask for several reasons.

AirFit P10 includes small, medium and large nasal pillows cushions, and it comes with a standard size headgear. It is also available in a for Her variant with smaller size option and uniquely styled headgear. Compatible with AirMini™.

- 1 **Adjustable QuickFit headgear** with the split-strap design is easy to put on, take off and adjust.
- 2 **Nasal pillows with dual-wall technology** seals comfortably on contact.
- 3 **QuietAir™ woven-mesh venting** diffuses air gently to reduce night-time disruptions.





## AirFit™ N30i

Nasal cradle mask

### Freedom

ResMed AirFit N30i features a top-of-the-head tube design that makes it easier to sleep in many positions and closer to your bed partner. It's a great option if you sleep on your side or stomach.

AirFit N30i nasal cradle cushions are interchangeable with AirFit P30i nasal pillows cushions. AirFit N30i is available in small and standard frame sizes. Both come with small, medium and small-wide cushions.

- 1 **SpringFit frame** adapts to different facial profiles to provide a personalised fit.
- 2 **Soft nasal cradle cushion** with QuietAir diffused venting to enhance comfort and provide a quieter night's sleep.
- 3 **No front-end tube** means the mask stays out of sight and out of the way.



## AirFit™ P30i

Nasal pillows mask


### Freedom

With a top-of-the-head tube design and a nasal pillows cushion, AirFit P30i makes it easier for you to sleep in many different positions and with less mask on your face.

AirFit P30i nasal pillows cushions are interchangeable with AirFit N30i nasal cradle cushions. AirFit P30i is available in small and standard frame sizes. The small starter pack comes with small and medium cushions; the standard starter pack comes small, medium and large size cushions.

- 1 **SpringFit frame** adapts to different facial profiles to provide a personalised fit.
- 2 **QuietAir vent** is designed to minimise air flow and night-time disruptions.
- 3 **No front-end tube** means the mask stays out of sight and out of the way.





## AirFit™ N20

Nasal mask

### Universal Fit


AirFit N20 is designed to fit a broad range of face shapes and sizes. Thanks to its under-eye frame for greater visual freedom and all-round comfort, AirFit N20 has become our most popular nasal mask since its release.

AirFit N20 is available in medium and large sizes. It is also available in a for Her variant with smaller size option and uniquely styled headgear. Compatible with AirMini.

- 1 InfinitySeal™ silicone cushion is designed to prevent blow-outs and seal well under variable therapy pressures.
- 2 Soft frame padding has been incorporated into the face-side of the frame to help provide comfort at every touchpoint.
- 3 Magnetic clips guide the headgear to the frame quickly and easily.\*



\*Use of masks with magnetic components is contraindicated in patients with the following pre-existing conditions: a metallic hemostatic clip implanted in the head to repair an aneurysm and/or metallic splinters in one or both eyes.



## AirTouch™ N20

Nasal mask

### Ultra Soft

AirTouch N20 features an UltraSoft™ memory foam cushion, which is designed to provide a light, breathable seal. It shares the same frame and headgear as AirFit N20 so users can swap between cushions if desired.

AirTouch N20 is available in medium and large sizes. It is also available in a for Her variant with smaller size option and uniquely styled headgear.

- 1 UltraSoft memory foam is a softer alternative than silicone and designed to reduce the risk of red marks and skin irritation.
- 2 A disposable cushion (to be replaced monthly) helps make maintaining hygiene and cleanliness easier than ever before.
- 3 Magnetic clips guide the headgear to the frame quickly and easily.\*



**14 DAY** REPLACEMENT GUARANTEE

If you're not completely satisfied with the UltraSoft memory foam cushion, you may exchange the cushions that are a part of the AirTouch N20 Starter Kit for an AirFit N20 InfinitySeal™ silicone cushion within 14 days of purchase. T&Cs apply. See in store for details.



# BECOME A VIP TODAY...

Join sleepvantage and be treated to a range of exclusive member benefits.\*

- ✓ **Extended warranty\*\***  
Additional sleepvantage membership warranty on ResMed products.
- ✓ **Feel supported**  
Access to our friendly product support line.
- ✓ **Stay informed**  
Access to a range of educational articles and videos.
- ✓ **Get special offers**  
Be part of exciting promotions and special offers.

For more information, contact sleepvantage today



[sleepvantage.com.au](https://sleepvantage.com.au)



1300 305 705 (AUS) or 0800 737 633 (NZ)



[info@sleepvantage.com.au](mailto:info@sleepvantage.com.au)

\*sleepvantage membership is available to residents of Aus and NZ.

\*\*Terms and conditions apply. See website for more information.



#### ResMed Ltd

1 Elizabeth Macarthur Drive  
Bella Vista NSW 2153 Australia  
1300 305 705 (AUS) or 0800 737 633 (NZ)

AirFit, AirMini, AirTouch, InfinitySeal, QuickFit, QuietAir, sleepvantage and UltraSoft are trademarks and/or registered trademarks of the ResMed family of companies. For patent information, see [ResMed.com/ip](https://ResMed.com/ip). Specifications may change without notice.

© 2021 ResMed Ltd. 10111537r4 2021-04