


# **Awaken**



#### Sleep therapy made easy

ResMed understands that caring for your mask and equipment can be time-consuming, confusing and can come with unexpected costs. Our range of therapy plans is designed to make payments more manageable and better match your lifestyle needs.

- ✓ The latest gear for optimal therapy delivered to your door
- ✓ Up to 5 years warranty on your device\*
- ✓ sleepvantage<sup>™</sup> membership and exclusive support<sup>\*</sup>
- ✓ Easy fixed weekly payments over 36 months
- AirMini App available for therapy tracking and coaching

#### About AirMini

AirMini is the world's smallest CPAP device that's packed with our latest technologies, making it suitable for regular therapy at home and while you're away. With the same established therapy algorithms as our best-selling bedside devices, you'll never have to choose between freedom and performance again.

- Three therapy modes: AutoSet<sup>™</sup>, AutoSet for Her and CPAP
- HumidX<sup>™</sup> waterless humidification for portable comfort\*\*
- AutoRamp<sup>™</sup> with sleep onset detection makes the beginning of each therapy session more comfortable
- Expiratory pressure relief (EPR™) reduces pressure during exhalation
- Dedicated AirMini App for easy setup and ongoing support
- Bluetooth<sup>™</sup> connectivity for remote therapy control from your smart device via AirMini App
- Designed and optimised for AirFit<sup>™</sup> 30, AirTouch<sup>™</sup> 20, AirFit 20 and AirFit P10 for AirMini masks.

### AirMini therapy plan options

Our AirMini plans are designed to make your therapy experience as seamless as possible. The table below describes what items you can expect to receive in each plan and how often you'll receive them.

	AirMini <b>Basic</b>	AirMini <b>Premium</b>	
AirMini device AirFit mask* + AirMini setup pack** AirMini filters AirMini mount system AirMini travel bag Mask headgear Mask cushion CPAP wipes tub + Plus one of the following: HumidX AirFit mask* + AirMini connector**	6	6	Starter pack
AirFit mask* AirMini filters CPAP wipes tub AirMini tubing	6 - 1	6 6 1	12 & 24 months
Mask headgear Mask cushion	-	1 1	18 & 30 months
*AirMini compatible and excludes AirTouch mask systems. **Not applicable for AirFit P10 and AirFit N30 for AirMini mask packs.	All this value for \$14/week (Total over 36 months \$2,283)	All this value for \$17/week (Total over 36 months \$2.751)	

<sup>\*</sup>sleepvantage membership required. Members receive an extra 3 years warranty in addition to the 2-year manufacturer's warranty

<sup>\*\*</sup>Currently available for all AirMini compatible masks excluding AirFit F30.

<sup>&</sup>lt;sup>†</sup>Based on volume comparisons done by ResMed internal testing prior to the time this piece was printed.

### ResMed AirMini therapy plans

Our therapy plans are designed to ensure you always have the latest products and support you need for ongoing therapy success – right from the very first night.





## AirMini App puts therapy in your hands

Simply download the AirMini App to your smart device, pair using Bluetooth wireless technology and enjoy convenient access to your comfort settings. AirMini App also provides you with a 30-day therapy history and therapy score each morning, so you can stay informed and motivated.

